




# Vision Monitor

## MonBaby *Box*



The **MonBaby** *Box* is a stand-alone system for recording flash ERGs and VEPs.



### Mechanical specifications

1/1

For one optical head:

Width	400 mm
Height	200 mm
Depth	370 mm
Weight	9.5 kg

### Electrical specifications

230 V, 0.7A or 110V, 1.4A 50 or 60 Hz

The **MonBaby** *Box* contains its own isolation transformer.

It includes one power-outlet that can be used to power a laptop.

It can also be connected to a standard PC tower with a color monitor but, in that case, it requires an external isolation transformer.



**WARNING:** Do not use it in the presence of flammable anesthetics.