



Endeavor™ CR – Proven Performer in the OR

Reliable, Feature Rich, Rugged, Portable

The Endeavor CR is the true multi-modality intraoperative monitoring system. In spinal, cranial, vascular or orthopedic procedures, the OR Staff needs to know that neurological and physiologic systems are not being compromised.

Feature Rich

- Capability to store virtually unlimited number of procedure templates
- Configurable to meet any IOM requirement by providing 4 independently controlled high-level stimulators, 12 stimulator sites and 1 dedicated low-level stimulator


Data Capture

- Critical data can be captured during procedure with Multiple Data Type software (MDT) that provides a comprehensive set of recording parameters available simultaneously
- Review data with different filters and montages with Natus Neurology's unique signal processing
- Ability to store all data during a case or only representative samples


Rugged & Portable

- Fits in the overhead storage bin on most commercial aircraft
- Able to withstand demanding OR environments





NicVue Patient Administration software



Acquisition Software


Connectivity – Remote Access

- NicVue Connect HL7 package seamlessly connects diagnostic equipment to the Hospital Information System (HIS)
- Bi-directional interface can help eliminate manual data entry, reduce transcription errors and save valuable time by importing patient demographics from your HIS
- NicVue Connect can assist with the centralization of clinical reports facilitating the move to a paperless healthcare model
- Remote access to support physician interaction using several supported communication packages

Supplies & Accessories

Personalize your Endeavor CR while promoting patient comfort with our broad portfolio of neurodiagnostic supplies and accessories. Healthcare professionals count on Natus Neurology to help them deliver exceptional cost-conscious care with reliable results.


Our dedicated customer service team provides a streamlined ordering and shipping process to save you time and money.



Subdermal Needle Electrodes



Disposable Corkscrew Electrode



Disposable Adhesive Electrodes



Service

Natus Neurology is committed to providing exemplary service to our customers. Our dedicated and experienced Customer Service Team will assist with every aspect of an order. To support our products, we provide factory-trained Field Technicians and Clinical Application Specialists for onsite support. Additionally, we provide an in-house Technical Support Team, staffed with experts, and a strong distribution network in International Markets to offer a wide range of service options. Allowing our customers more time to care for their patients is our goal. Customer loyalty is our reward.

Supplies

Natus Neurology offers a full range of neurodiagnostic accessories and supplies promoting patient comfort. Our dedicated customer service team provides a streamlined order and shipping process to save you time and money.

To learn more about Natus Neurology Service Programs or our full line of Supplies and Accessories, contact your local distributor or sales representative.
US Customers Call: **1-888-662-7632**



Natus Neurology Incorporated

3150 Pleasant View Road
Middleton, WI 53562 USA
Tel: 1-800-356-0007
1-608-829-8500
Fax: 1-608-829-8709
www.natus.com

© 2014 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.